EasyTom CT Scanner 3D Micro Computed Tomography & Digital Radioscopy System

CT Scanner with Ultra 3D Accuracy and Resolution for Large Volume Inspection

EasyTom 130 & 150 Micro CT Scanners are x-ray inspection machines allowing a collection of complex internal and external geometry. The EasyTom 3D CT scanners features high-resolution digital radioscopy, versatility for a wide variety of applications, 6 motion axis and large volume inspection.



Attain measurements in real time with high resolution digital radioscopy

Attains Highly Precise Measurements

- Capture highly accurate measurements at +/- $10\mu m$ Accuracy with Resolution to $4\mu m$
- High accuracy motorized rotation and 3 Axis translations


Provides Easy-to-Use 3D Scanning Capabilities

- Easy to use 3D micro and computed tomography system

Versatile for a Variety of Applications and Parts

- Easily verify structure of 3D printed metal parts
- Attain external and Internal Surface Geometry output as .STL for use with popular 3D Scan Data Processing software
- Add an optional 3D Visualization and post processing software suite to fit any application: Inspection, Reverse Engineering, Analysis, Porosity, Fiber Alignment, Wall Thickness, Comparison to CAD 3D Color Maps and much more
- Full inspection of large size samples and high volume inspection.
- Full inspection volume (ØxH)*: 32 cm (12.6") x 42 cm (16.5")





Part Capability

- Full inspection volume 32 cm (12.6") diam* 42 cm (16.5")
- W = 410 mm, using the system in 3D CT and 2D Radio: Max. Shift 2: 460 mm max
- $H = 410 \, \text{mm}$
- SDD (distance between the emitter and the detector) = 900 mm (max. available space)



Safety Cabinet	 Footprint: 2100x1100x2000 mm / 82.6"x43.3"x78.7" Lead / Steel construction and X-ray safety interlocks, designed to meet X-ray safety regulations Motorized door with automatic locker during X-Ray emission. Large scanning volume (diameter x height): 320 mm x 420 mm
Mechanics	 High accuracy motorized rotation and translation axis. Imager lateral and vertical shift option for enlarged field of view and decreased ring artifacts. Air-bearing rotation stage option, takes up of sample weight.
X-Ray Generator	Several options and combinations available: • Sealed or open type micro-focus tube. • Open nano-focus tube (160 kV) • Voltage up to 230 kV (several options available). • Resolution down to 400nm/voxel • Various targets and filament types available.
lmager	Several options and combinations available: • High resolution flat panel detector • Large area flat panel detector • CCD sensor

RX SOLUTIONS CT SOFTWARE: X- ACT

CT Acquisition

CT Acquisition Modes: conventional, helical, stack, laminography, continuous or step by step rotation. Ergonomy: wizard mode for non experts, automation mode for single click acquisition to inspection workflow. Radiography filter enhancement, 2D video sequence acquisition, 3D measurements. Automatic black and gain calibration & sample repositioning.

CT Reconstruction

Real time artifacts corrections: focal spot drift, ring artifacts, beam hardening, phase contrast. Geometry compensation: automatic correction of the rotation center and other geometric parameters. Easy and intuitive 3D optimization of the reconstruction volume using test slices. On the fly reconstruction of a running acquisition.

Workstations

System-integrated acquisition workstation. Standalone reconstruction workstation with powerful GPU.

Analysis Software (Optional)

VGStudio or VGStudio MAX: 3D Visualization and post-processing software with metrology, CAD comparison, porosity, and wall thickness analysis module

**X-ACT Software is from RX Solutions SAS - VGStudio and VGStudio MAX is from Volume Graphics, Inc.

Manufactured by RX Solutions SAS, Chavanod, France





